



Non-Cuffed Wire Reinforced Endotracheal Tubes - Silicone

The transparent silicone tube allows easy visual inspection. The metal spiral wire enables the position of the tube to be easily confirmed through x-ray control. This special design prevents the tube from kinking whilst it is in use. The bevel tip allows easy insertion. Intended for VETERINARY USE ONLY.



Order No.	Tube I.D. (mm)	Tube O.D. (mm)	French Size	Length (cm)
V-PATR-25	2.5	4.0	12	14
V-PATR-30	3.0	4.7	14	16
V-PATR-35	3.5	5.3	16	18
V-PATR-40	4.0	6.0	18	20
V-PATR-45	4.5	6.7	20	22
V-PATR-50	5.0	7.3	22	24

Cuffed Wire Reinforced Endotracheal Tubes - Silicone

The transparent silicone tube allows easy visual inspection. The metal spiral wire enables the position of the tube to be easily confirmed through x-ray control. This special design prevents the tube from kinking whilst it is in use. Low pressure cuff enables atraumatic positioning. The bevel tip allows easy insertion. Intended for VETERINARY USE ONLY.

Order No.	Tube I.D. (mm)	Tube O.D. (mm)	French Size (mm)	Length (cm)
V-PETR-6	6.0	8.7	26	28
V-PETR-7	7.0	10.0	30	30
V-PETR-8	8.0	11.3	34	32
V-PETR-9	9.0	12.7	38	34
V-PETR-10	10.0	14.0	42	36

