

ACCESSORIES

Stopcocks



Order No.

V-POWSC-FLL-MLL-LP

Remarks

One-way plastic stopcock, female Luer lock to male Luer lock.



V-PTWSC-2FLL-MLL-R-LP

Three-way plastic stopcock, two female Luer locks to male Luer lock with right side arm.

Multi-Purpose Tubing Adaptors



Order No.

V-PFLLA-VTA


Remarks

Female Luer lock to tapered catheter end. Fitting accepts 10-20 French inside diameter.



V-PFLLA-VTA-L

Female Luer lock to large tapered catheter end. Fitting accepts 20-40 French inside diameter.



V-PMLLA-VTA-L

Male Luer lock to large tapered catheter end. Fitting accepts 20-40 French inside diameter.



V-FD-100

Silicone fixation disc used to stabilise percutaneous drainage catheters.



V-P79

7.0-9.0mm double-ended tapered connector.

Multi-Purpose Tubing Adaptors

Order No.

V-P412

Remarks

4.0-12mm double-ended tapered connector.



Injection Ports

Order No.

V-VIP-200

Remarks

Luer lock cap to allow injection of solutions through cap



Clamps

Order No.

V-CLAMP-W

Remarks

Clamp to occlude drainage catheter tubing.

V-CLAMP-W-L

Large clamp to occlude drainage catheter tubing



Silicone Adhesive

Order No.

V-ADH-2

Remarks

Single tube of silicone adhesive. Recommended for repair of balloon cuffs



Luer Lock Plug

Order No.

V-PFLLP-100

Remarks

Luer lock plug to occlude female Luer lock.

