




Manifolds without Sentinel DO NOT meet the requirements of C.S.A. Porter Instrument Co., Inc. accepts no responsibility for the installation and/or use of manifolds without Sentinel.




5201K Desk Sentinel™ Add-On Kit for 3022MX Manifold
 includes:
 5201 Desk Sentinel
 C-394-000 (2) Hi/Lo Pressure Switch
 A-2041 Sentinel Interconnect Board Assembly




5251K Wall Sentinel™ Add-On Kit for 3022MX Manifold
 includes:
 5251 Wall Sentinel
 C-394-000 (2) Hi/Lo Pressure Switch
 A-2041 Sentinel Interconnect Board Assembly



7000 Oxygen Regulator, Hi Pressure Gauge Only
7500 Nitrous Oxide Regulator, Hi Pressure Gauge Only




***A-2144 diss/Copper Pigtails (straight)**
 Oxygen-1/2" O.D. copper x diss O₂ Female,
 Nitrous Oxide-3/8" O.D. copper x
 diss N₂O Female



***A-2145 diss/Copper Pigtails (90° elbow)**
 Oxygen-1/2" O.D. copper x diss O₂ Female,
 Nitrous Oxide-3/8" O.D. copper x
 diss N₂O Female

***Use for pre-installation of Manifold when piping is concealed in the tank room (see template FM-188)**



3525 Set of two (2) O₂/N₂O Tank Restraints

PA-514-000 **O₂ Pipeline Label, vinyl**
 PA-515-000 **N₂O Pipeline Label, vinyl**