


# CHEST DRAINAGE TUBES – SILICONE



## DESCRIPTION

- Open end single lumen catheter with large multiple sideholes over the distal portion of the catheter.
- The catheter can be used for the drainage of air or fluid from body cavities including the pleural space of the thoracic cage.
- The catheter is placed by blunt dissection. Metal stylets available to stiffen catheter and assist placement.
- Markings at 2.0 and 4.0 cm from the proximal sidehole assist in correct catheter placement.
- Intended for VETERINARY USE ONLY.
- **Maximum sterilisation temperature 134°C for 3 minutes.**

## Order No.

## Remarks

V-PCT-20	20 French straight drainage tube 42 cm long with 4.0 mm I.D. and 6.6 mm O.D. Sideholes over distal 10 cm of catheter.
V-PCT-30	30 French straight drainage tube 42 cm long with 7.0mm I.D. and 10 mm O.D. Sideholes over distal 12 cm of catheter.

## METAL INTRODUCTION STYLET FOR USE WITH CHEST DRAINAGE TUBES



## Order No.

## Remarks

V-TCS-119-45	Metal stiffening stylet, 3.2 mm O.D., 45 cm (18") long. Recommended for use with 20 French Chest Drainage Tubes.
V-TCS-187-45	Metal stiffening stylet, 4.75 mm O.D., 45 cm (18") long. Recommended for use with 30 French Chest Drainage Tubes.

