




# THAL-QUICK CHEST DRAINAGE CATHETER SETS

## WITH SELDINGER INTRODUCTION TECHNIQUE



### DESCRIPTION

- Used for pleural fluid drainage.
- Chest tube placement is facilitated utilising the Seldinger percutaneous entry technique. Refer to page 13 for instructions.
- The chest tube inserter, fitted to the chest tube, provides a smooth transition for ease of insertion over the wire guide.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

### Order No.

### Remarks

V-TQTS-1400	14 French, 41 cm (16") long. Four large sideports are positioned at distal end of chest tube.
V-TQTS-1800	18 French, 41 cm (16") long. Four large sideports are positioned at distal end of chest tube.
V-TQTS-2400	24 French, 41 cm (16") long. Four large sideports are positioned at distal end of chest tube.
V-TQTS-2800	28 French, 41 cm (16") long. Four large sideports are positioned at distal end of chest tube.

### PLEASE NOTE

*V-PFLA-VTA-L may be used with this product to adapt to a Luer lock fitting, refer to page 34. One-way Heimlich valve, order no. V-CASP-A-FORD may also be used in conjunction with this product, refer to page 14.*