




HEIMLICH VALVE



HEIMLICH CHEST DRAINAGE VALVE



DRAINAGE CHANNELLING TUBE


ORDER No.

Remarks

V-CASP-A-FORD
V-CTU14.0-30

Heimlich one-way valve and connecting tube.
Connecting tube only.

ABDOMINOCENTESIS CATHETER



CATHETER



OBTURATOR

DESCRIPTION

- Abdominocentesis catheter set containing a single lumen open ended catheter with distal sideholes and proximal Luer lock fitting.
- Blunt obturator to assist the safe placement of the catheter.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

Order No.

Remarks

V-EACC-10.2-10-BO

Set contains a 10.2 French, 10 cm (4") long catheter with sideholes and blunt obturator.