


# LEPAGE BALLOON DRAINAGE CATHETER



## DESCRIPTION


- A closed ended single lumen catheter with sideholes, a proximal balloon and a moveable silicone retention disc.
- The catheters may be used for the drainage of the thoracic and abdominal cavities.
- Supplied in peel-open packages.
- Intended for VETERINARY USE ONLY.
- **Maximum sterilisation temperature 115°C for 20 minutes.**

Order No.	French	Length	Remarks
V-EBDC-20-30.0-LEPAGE	20	30 cm	Single lumen catheter with sideholes.

## REFERENCES

*Dr Olivier Lepage, Earle Nationale Veterinaire de Lyon, France.*

## METAL INTRODUCTION STYLET FOR USE WITH LEPAGE BALLOON DRAINAGE CATHETER



STIFFENING STYLET  
Stainless Steel

Order No.	Remarks
V-TCS-90-40	Metal stiffening stylet, 2.3 mm O.D., 40 cm long. Recommended for use with 20 French drainage catheters.