




'FLEXI-DRAIN' MULTI-PURPOSE WOUND DRAIN – SILICONE



DESCRIPTION

Order No.

V-PWD-12-40



Top View



End View

- A silicone tubular drain that may be trimmed or used as one piece for surgical wound drainage. The drain can be split along its length or trimmed to length to suit individual requirements.
- Drain material is non-radiopaque.
- Intended for VETERINARY USE ONLY.
- *Maximum sterilisation temperature 134°C for 3 minutes.*

Remarks

40 cm (16") long drain consisting of twelve catheter lumens joined together, each of 3.0 mm diameter.