



MULTIPLE LUMEN JUGULAR CATHETER SETS

PEEL-AWAY® SHEATH INTRODUCTION TECHNIQUE



DESCRIPTION

- Used for venous pressure monitoring, blood sampling and administration of drugs and fluids.
- The catheters incorporate separate, non-communicating vascular access lumens.
- The multiple lumen design reduces the necessity of multiple venipunctures.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

ORDER No. Remarks

DOUBLE LUMEN CATHETER SET

V-PPC-UDLM-751-LSC 7.5 French double lumen catheter, 25 cm (10") long with a 13 gauge, 7.0 cm long Peel-Away® sheath needle and a 10 ml syringe.

TRIPLE LUMEN CATHETER SET

V-PPC-UTLM-701-LSC 7 French triple lumen catheter 25 cm (10") long with a 13 gauge, 7.0 cm long Peel-Away® sheath needle and a 10 ml syringe.