


MULTIPLE LUMEN JUGULAR CATHETER SETS

SELDINGER INTRODUCTION TECHNIQUE



DESCRIPTION

- Used for venous pressure monitoring, blood sampling and administration of drugs and fluids.
- The catheters incorporate separate non-communicating vascular access lumens.
- The multiple lumen design reduces the necessity of multiple venipunctures.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

ORDER No. Remarks

DOUBLE LUMEN CATHETER SET

V-UDLM-751J 7.5 French, 13 gauge catheter, 15 cm (6") long with .035" (0.89 mm) straight and 3.0 mm 'J' tip wire guide.

V-UDLM-751J-LSC 7.5 French, 13 gauge catheter, 25 cm (10") long with .035" (0.89 mm) straight and 3.0 mm 'J' tip wire guide.

TRIPLE LUMEN CATHETER SET

V-UTLM-701J-LSC 7 French, 13 gauge catheter, 25 cm (10") long with .035" (0.89 mm) straight and 3.0 mm 'J' tip wire guide.