




# DISPOSABLE CYTOLOGY BRUSH



## DESCRIPTION

- Used to collect cells from the bronchi or gastrointestinal tract for cytology.
- Ring handle enables the brush to be retracted into the sheath to prevent contamination of the sample.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

## ORDER No.

## Remarks

V-ECB-5-180-3-S

One 5 French, 180 cm long cytology brush with 3.0 mm brush diameter.