



## VET BIOSIST™ OCULAR DISC

### DESCRIPTION

- Applications: *VET BIOSIST™* is for use as a surgical patch for soft tissue repair or reinforcement. It promotes wound healing and tissue remodelling by providing a biocompatible, absorbable scaffold for tissue ingrowth.
- Suitable applications are corneal lacerations, periodontal grafts, delayed closure of dermal wounds, organ patching and adhesion barrier.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

### ORDER No.

### Remarks

V-LSIS-D-010-4

Four ocular discs 10 mm in diameter, individually packed.



## VET BIOSIST™ MULTILAMINATE SHEET

### DESCRIPTION

- Applications: *VET BIOSIST™* is for use as a surgical patch for soft tissue repair or reinforcement.
- Multilaminated sheets are suitable for applications where high physical pressures or chemical effects of body fluids demand a stronger patch than the standard single layer such as thoracic cage, muscular defects, and large animal skin lesions.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

### ORDER No.

### Remarks

V-ASIS-M4-070100

One multilaminated *VET BIOSIST™* sheet 70 mm x 100 mm.