



## STOPCOCKS

ORDER No.

Remarks

V-POWSC-FLL-MLL-LP

One-way plastic stopcock, female Luer lock to male Luer lock.



ORDER No.

Remarks

V-PTWSC-2FLL-MLL-R-LP

Three-way plastic stopcock, two female Luer locks to male Luer lock with right side arm.



ORDER No.

Remarks

V-P2V

One-way plastic stopcock, female Luer slip to male Luer slip.



ORDER No.

Remarks

V-P3V

Three-way plastic stopcock, two female Luer slips to one male Luer slip.




## MULTI-PURPOSE TUBING ADAPTORS

ORDER No.

Remarks

V-PFLLA-VTA

Female Luer lock to tapered catheter end; fitting accepts 10-20 French inside diameter.



ORDER No.

Remarks

V-PFLLA-VTA-L


Female Luer lock to large tapered catheter end; fitting accepts 20-40 French inside diameter.

# ACCESSORIES



## MULTI-PURPOSE TUBING ADAPTORS


ORDER No.	Remarks
V-PMLLA-VTA-L	Male Luer lock to large tapered catheter end; fitting accepts 20-40 French inside diameter.



ORDER No.	Remarks
V-P412	4.0-12 mm double-ended tapered connector.



ORDER No.	Remarks
V-P79	7.0-9.0 mm double-ended tapered connector.



## FIXATION DISC

ORDER No.	Remarks
V-FD-100	Silicone fixation disc used to stabilise percutaneous drainage catheters.



## INJECTION PORTS

ORDER No.	Remarks
V-VIP-200	Luer lock cap to allow injection of solutions through cap.



ORDER No.	Remarks
V-VIP-300	T-piece injection port. Allows connection between Luer lock fittings.



## INJECTION PORTS

ORDER No.	Remarks
-----------	---------

V-BDV-100	Bi-directional valve.
-----------	-----------------------



## LUER LOCK PLUG

ORDER No.	Remarks
-----------	---------

V-PFLLP-100	Luer lock plug to occlude female Luer lock.
-------------	---------------------------------------------



## CLAMPS

ORDER No.	Remarks
-----------	---------

V-CLAMP-W	Clamp to occlude drainage catheter tubing.
-----------	--------------------------------------------



ORDER No.	Remarks
-----------	---------

V-CLAMP-W-L	Large clamp to occlude drainage catheter tubing.
-------------	--------------------------------------------------



## SILICONE ADHESIVE

ORDER No.	Remarks
-----------	---------

V-ADH-2	Single tube of silicone adhesive for use when repairing balloon cuffs.
---------	------------------------------------------------------------------------

# CONVERSION CHART

## METRIC/IMPERIAL MEASUREMENT


### LENGTH

METRIC		IMPERIAL
1 millimetre (mm)		0.0394 inches (in)
1 centimetre (cm)	10 mm	0.3937 inches (in)
1 metre (m)	100 cm	1.0936 yards (yd)


  

IMPERIAL	METRIC	
1 inch (in)	2.54 centimetres (cm)	
1 foot (ft)	12 in	0.3048 metres (m)
1 yard (yd)	3 ft	0.9144 metres (m)

### TEMPERATURE



### FRENCH CATHETER SCALE



### STUBS NEEDLE GAUGE

