


DOUBLE SPIKE WIDE BORE INFUSION SET



DESCRIPTION



- A complete double spike tubing set used mainly for adult horses in stalls.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

Order No.

Remarks

V-LACT-24-450-2S-C-EC4 Double wide bore spike, 24 French (8.0 mm) 450 cm total tubing length, large in-line chamber, roller clamp, injection port and rotating cuff male Luer lock.

SINGLE SPIKE WIDE BORE INFUSION SET



DESCRIPTION

- A complete single spike tubing set used mainly for adult horses in stalls.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

Order No.

Remarks

V-LACT-24-450-S-C-EC2 Single wide bore spike, 24 French (8.0 mm) 450 cm total tubing length, large in-line chamber, roller clamp, injection port and rotating cuff male Luer lock.