


WIDE BORE INFUSION SET (NO CHAMBER)



DESCRIPTION


- A simple single spike set for rapid fluid administration.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

Order No.

Remarks

V-LACT-24-200-S-NC-CC1 Single wide bore spike, 24 French (8.0 mm) 200 cm length tubing with rotating cuff male Luer lock.

BAG TO BAG TRANSFER SET



DESCRIPTION

- A simple double spiked connector to allow additional fluid bags to be attached to an initial single spike set.
- Supplied sterile in peel open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

Order No.

Remarks

V-LACT-24-40-2S-NC-EC3 Two wide bore spiked accessory tubing, 24 French (8.0 mm) 40 cm length tubing.