



FOLEY CATHETERS – SILICONE




DESCRIPTION

- A closed endhole catheter with distal balloon cuff used for intermittent drainage of urine from the bladder.
- The balloon aids in retaining the catheter inside the bladder.
- Intended for VETERINARY USE ONLY.
- **Maximum sterilisation temperature 115°C for 20 minutes.**

Order No.	Remarks
V-PFC8-30	8 French Foley, 30 cm (12") long with 1.3 mm I.D. and 2.7 mm O.D.
V-PFC8-55	8 French Foley, 55 cm (22") long with 1.3 mm I.D. and 2.7 mm O.D.
V-PFC10-30	10 French Foley, 30 cm (12") long with 2.0 mm I.D. and 3.3 mm O.D.A
V-PFC10-55	10 French Foley, 55 cm (22") long with 2.0 mm I.D. and 3.3 mm O.D.
V-PFC12-30	12 French Foley, 30 cm (12") long with 2.2 mm I.D. and 4.0 mm O.D.
V-PFC12-55	12 French Foley, 55 cm (22") long with 2.2 mm I.D. and 4.0 mm O.D.
V-PFC14-55	14 French Foley, 55 cm (22") long with 2.5 mm I.D. and 4.7 mm O.D.
V-PFC16-55	16 French Foley, 55 cm (22") long with 2.8 mm I.D. and 5.3 mm O.D.
V-PFC18-55	18 French Foley, 55 cm (22") long with 3.2 mm I.D. and 6.0 mm O.D.
V-PFC20-65	20 French Foley, 65 cm (26") long with 3.6 mm I.D. and 6.7 mm O.D.
V-PFC22-65	22 French Foley, 65 cm (26") long with 4.0 mm I.D. and 7.3 mm O.D.
V-PSFC28-75	28 French Foley, 75 cm (30") long with 6.0 mm I.D. and 9.3 mm O.D.
V-PSFC28-100	28 French Foley, 100 cm (40") long with 6.0 mm I.D. and 9.3 mm O.D.



Inflate balloon with sterile water only.



Catheters of size 14 gauge and upwards may be placed using a rigid stylet.



Luer lock adaptor V-PFLA-VTA-L may be used with the Foley catheter, refer to page 34.