


# GUTTURAL POUCH ENDOSCOPIC FLUSHING CATHETER



## DESCRIPTION



- An endoscopically guided flushing catheter for lung lavage or introduction into the guttural pouch to flush out pouch debris.
- Visual indicators at 2.5, 5.0 and 7.5 cm from the distal tip, aid in introduction.
- Supplied sterile in peel open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.


Order No.	Remarks
-----------	---------

V-GPFC-6.0-NT-250-P-NS-O	6 French, 250 cm long PTFE catheter with distal tip markings.
--------------------------	---

V-GPFC-7.0-NT-250-P-NS-O	7 French, 250 cm long PTFE catheter with distal tip markings.
--------------------------	---



# UTERINE FLUSHING CATHETER – SILICONE



## DESCRIPTION

- A single lumen cuffed catheter for the administration of therapeutic flushing solutions into the uterus.
- The balloon cuff can be inflated to seal the internal os of the uterine cavity and prevent fluid loss.
- The proximal dual fitting allows attachment for both irrigation and drainage.
- Repair balloon cuff and inflation line are available.
- Intended for VETERINARY USE ONLY.
- **Maximum sterilisation temperature 115°C for 20 minutes.**

Order No.	Remarks
-----------	---------

V-PUF8-80	Uterine flushing tube, 33 French, 80 cm (32") long with 8.0 mm I.D. and 11 mm O.D.
-----------	--

V-CLAMP-W-L	Large clamp for the occlusion of drainage catheters and tubing replacements for the above set. This item is not autoclavable. (Illustrated page 36).
-------------	--